

## RTI-2101 Simple Test Station

The RTI –2101 Simple Test Station is a small, yet reliable and stable test station.


### Specifications

**4" Stainless Steel Chuck**  
**Vacuum Rings x 1 for wafer 2" and single die**  
**Side Driven Chuck Stage X-Y 3"-3"**

**Resolution 10 micron**  
**Platen Bar x 2 (East/West)**  
**Chuck Theta 30 degrees**

- Capable for DC, RF, HVHC
- Portable and can be moved between laboratories
  - o Share instruments if they're not accessible in your lab!
- Placed on top of an optical table
  - o Use directly under a solar simulator
- Small test station
  - o For chips or dies
  - o One-off testing away from production testing
  - o Classroom use for teaching and tutorial

Digital Microscope Quoted Upon Request



Optics quoted upon request

# RTI-2101 Simple Test Station

## **Probe Positioners –**

- Linear X-Y-Z Travel 12mm-12mm-12mm
- 40 Thread/Inch
- Stainless Steel Lead Screw
- Magnet Base
- Right and Left Hand Use

## **Tip Holders -**

- Spring Type Tip Holder
- Coaxial cable 1.5 m long terminated with BNC male

## **Tungsten Tip –**

- 5 micron dia.
- Box of 10

---

---

## **Dimensions -**

320mmWx320mmDx320mmH

---

---

For ordering information contact:  
Michelle Bell  
Radiant Technologies, Inc  
2835 Ste B Pan American Fwy NE  
Albuquerque, NM 87107  
[www.ferrodevices.com](http://www.ferrodevices.com)  
Tel - 505-842-8007  
Fax - 505-842-0366